


Active Vibration Isolation Platform SW1100x900-4.0



- dimensions: 1100x900x130mm
- weight: 144kg
- maximum load capacity: 600kg
- material of housing: aluminum
- surface: anodized
- configuration: 4 internal isolation elements (16 sensor/actuator pairs)
- electronics: 4-channel controller
- cables: 4 D-sub cables (3m long)

D-sub connectors require a space of 100mm



position of D-sub connectors

4-channel controller



		Allgemein-toleranzen ISO 2768-f		Werkstück-kanten DIN ISO 13715		Werkstoff:		Maßstab:	
		Datum		Name		Active Vibration Isolation Platform SW1100x900-4.0 and Controller			
Gezeichnet		07.08.2017		U. Stöbener					
Kontrolliert									
Norm						Zeichnungs-Nr.:			
Artikel-Nr.:									
Status	Änderungen	Datum	Name	F:\CAD\Technische Zeichnungen\MOD-Sandwich Zeichnungen\Sondermaß 1100x900\SW1100x900_4_Controller.idw					

A3