



**Active Vibration Isolation Platform SW1100x900-6.0**

- dimensions: 1100x900x130mm
- weight: 157kg
- maximum load capacity: 900kg
- material of housing: aluminum
- surface: anodized
- configuration: 6 internal isolation elements (24 sensor/actuator pairs)
- electronics: 2-channel controller & 4-channel controller
- cables: 6 D-sub cables (3m long)



		Allgemein-toleranzen ISO 2768-f		Werkstück-kanten DIN ISO 13715		Werkstoff:		Maßstab:	
		Datum		Name		<b>Active Vibration Isolation Platform SW1100x900-6.0 and Controller</b>			
Gezeichnet		07.08.2017		U. Stöbener					
Kontrolliert									
Norm						Artikel-Nr.:		Zeichnungs-Nr.:	
Status		Änderungen		Datum		Name		F:\CAD\Technische Zeichnungen\MOD-Sandwich Zeichnungen\Sondermaß 1100x900\SW1100x900_6_Controller.idw	

A3