



Active Vibration Isolation Platform SW1100x900-8.0

- dimensions: 1100x900x130mm
- weight: 170kg
- maximum load capacity: 1200kg
- material of housing: aluminum
- surface: anodized
- configuration: 8 internal isolation elements (32 sensor/actuator pairs)
- electronics: 2 x 4-channel controller
- cables: 8 D-sub cables (3m long)

4-channel controller



		Allgemein-toleranzen ISO 2768-f		Werkstück-kanten DIN ISO 13715		Werkstoff:		Maßstab:	
		Datum		Name		Active Vibration Isolation Platform SW1100x900-8.0 and Controller			
Gezeichnet		07.08.2017		U. Stöbener					
Kontrolliert									
Norm						Artikel-Nr.:		Zeichnungs-Nr.:	
Status		Änderungen		Datum		Name		F:\CAD\Technische Zeichnungen\MOD-Sandwich Zeichnungen\Sondermaß 1100x900\SW1100x900_8_Controller.idw	

A3